

### Capsule, Capsulotomy, Vannas



#### G-19776 Capsulotomy

- Sharp tips
- Angled blades
- Cut length 7.5mm
- 9.0mm blades
- 2.0mm pivot width
- 80.0mm overall length
- Chromium plated stainless steel, reusable
- G-19777 cut length 8.5mm



#### G-19620 Tuebingen Model Capsulotomy

- Sharp tips
- Straight blades
- Cut length 7.0mm
- 9.0mm blades
- 1.8mm pivot width
- 88.0mm overall length
- Chromium plated stainless steel, reusable



#### G-19760 Vannas Capsulotomy

- Sharp tips
- Extremely delicate
- Straight blades
- Cut length 4.5mm
- 9.0mm blades
- 1.8mm pivot width
- 88.0mm overall length
- Chromium plated stainless steel, reusable



#### G-19765 Vannas Capsulotomy

- Sharp tips
- Extremely delicate
- Curved blades
- Cut length 4.5mm
- 1.7mm pivot width
- 80.0mm overall length
- Chromium plated stainless steel, reusable




#### G-19775 Vannas Capsulotomy

- Sharp tips
- Extremely delicate
- Angled blades
- Cut length 5.0mm
- 7.0mm blades
- 80.0mm overall length
- Chromium plated stainless steel, reusable



#### Benefits of GEUDER<sup>®</sup> Chromium Plated Surgical Instruments

- ✓ Surface Hardness
- ✓ Corrosion Protection
- ✓ Ease of Cleaning and Sterilization
- ✓ Reduced Friction
- ✓ Long Lasting
- ✓ 67+ Years of Experience




### Corneal and Corneal Section



#### G-19695 Castroviejo Corneal Left and Right

- Blunt tips
- Left blades, curved
- Cut length 7.5mm
- Lower blade 9.5mm
- Upper blade 9.0mm
- 105.0mm overall length
- Chromium plated stainless steel, reusable
- G-19690 right blades



#### G-19605 Corneal Universal

- Blunt tips
- Lightly curved blades
- Cut length 6.0mm
- 8.0mm blades
- 105.0mm overall length
- Chromium plated stainless steel, reusable
- G-19600 sharp tips



#### G-19715 Troutman Corneal Left and Right

- Blunt, delicate tips
- Left blades, curved
- Cut length 6.0mm
- 7.0mm blades
- Tip to pivot 9.5mm
- 100.0mm overall length
- Chromium plated stainless steel, reusable
- G-19710 right blades



Many more styles available. Most popular products are in stock.

Need something you don't see here?  
Please contact us for assistance.

Our expert Customer Resource team will provide you with fast, friendly, personalized service!



## Enucleation

### G-19520 Enucleation

- Blunt tips
- Curved blades
- Strong curve
- Tip to pivot 40.0mm
- 135.0mm overall length
- Chromium plated stainless steel, reusable
- G-19525 medium curve
- G-19530 light curve



## Iris


### G-19680 Barraquer Iris

- Blunt tips
- Angled blades
- 7.0mm blades
- Cut length 5.5mm
- 56.0mm overall length
- Chromium plated stainless steel, reusable



## IOL Removal/Explantation

### G-32996 Snyder-Osher IOL Scissors



- For cutting acrylic and silicone foldable IOLs, lower notched blade to create IOL traction
- 18 gauge / 1.2mm shaft, overall length 150.0mm, chromium plated stainless steel, reusable
- Long model, 4.5mm blades
- G-32998 short model, 3.5mm blades
- Use with G-32997 Snyder-Osher Forceps

## Tenotomy

### G-19485 Stevens Tenotomy

- Blunt tips
- Curved blades
- Cut length 19.0mm
- Tip to pivot 28.0mm
- 110.0mm overall length
- Chromium plated stainless steel, reusable
- G-19480 blunt tips, straight blades, cut length 18.0mm



## Eye/Drape

### G-19375 Curved Eye

- Sharp tips
- Curved blades
- Cut length 18.0mm
- Various other tip lengths (see Geuder Catalog)
- Chromium plated stainless steel, reusable
- G-19355 sharp tips, straight blades, cut length 17.0mm



## Strabismus

### G-19465 Strabismus

- Blunt tips
- Curved blades
- Short blade
- Tip to pivot 28.0mm
- Cut length 19.0mm
- 110.0mm overall length
- Chromium plated stainless steel, reusable
- G-19460 blunt tips with straight blades



## Utility and Westcott

### G-19751 Westcott Type Stitch

- Sharp tips
- Curved blades
- Cut length 14.0mm
- 120.0mm overall length
- Chromium plated stainless steel, reusable
- G-19750 blunt tips, curved blades



## Contact Us

**One Source**  
SERVICE AND SUPPORT

Personalized customer service  
In-house experts  
Best-of-best products  
Sales, repair, refurbish, replace  
Custom design & development  
Committed to quality patient care

MR. FEMTO  
MR. PHACO

**DK**  
Duckworth & Kent  
Precision Ophthalmic Instrumentation

**Geuder**  
Precision made in Germany

**MANI**

**STERIMEDIX**

**AMCON**

**ROUND'S**

**NWBOC**  
NWSB Certified Women Owned Business

**Crestpoint Management**  
**Crestpoint Ophthalmics**  
9100 Watson Road  
St. Louis, MO 63126 USA

GSA/FSS Contract  
V797P-4438B  
Ophthalmic Instruments  
SIN A8a & A8b  
Effective  
5/1/2011 to 4/30/2021

Published by Crestpoint Management. Manufacturer trademarks, part data, and copyrighted materials are property of their respective owners and are used with permission. Geuder product ranges in this publication are exclusive in the USA / NA.

Phone: (800) 756-7773  
Phone: (314) 849-7773  
Fax: (314) 849-7776  
E-mail: [custsvc@crestpointmgt.com](mailto:custsvc@crestpointmgt.com)  
Web: [www.crestpointmgt.com](http://www.crestpointmgt.com)

Copyright 2018 Crestpoint Management / Crestpoint Ophthalmics. All rights reserved. All information in this publication is subject to change without notice.

ZCC425 (10/2018)

**Crestpoint**  
OPHTHALMICS